



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CALIBRATOR LYOPHIL. IN URINE FOR TRICHLOROACETIC ACID

Code Z17016
TO USE WITH THE KIT CODE GC17010

PURPOSE

These urinary calibrators are used to calibrate the GC system for the quantitative measurement of analytes contained therein. Each analyte has a concentration that allows the construction of a curve on a point. The calibrators are freeze-dried human matrix and should be handled as if they were a real patient sample.

RECONSTITUTION


Remove the metal seal and rubber stopper from the vial. Add exactly 1 ml of HPLC grade H₂O in the vial. Replace the rubber stopper, shake and let stand for 5 to 10 minutes. Before use, mix by inverting the vial to dissolve the material until a homogeneous clear solution.

STORAGE AND STABILITY

The calibrators are stable for 36 months from the date of preparation if stored at 2-8 ° C. After reconstitution are stable 7 days at 2-8 ° C and 1 month at -20 ° C. Do not use after the expiry date.

PRECAUTIONS

These calibrators in human matrix should be treated with care and treated as potentially infectious.

LOT		

PACKAGING AVAILABLE:

- Z17016 CALIBRATOR IN URINE FOR TRICHLOROACETIC ACID

4 x 1 ml




N° 003

06/2024

CONCENTRATIONS

CODE

Z17016

<i>ANALITE</i>	<i>UNITS OF MEASUREMENT</i>	<i>VALUE</i>
TRICHLOROACETIC ACID	mg/l	

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR *IN VITRO* DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com

N° 003	06/2024
---------------	----------------