



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CALIBRATOR LYOPHIL. IN URINE FOR trans 1,2-CYCLOHEXANEDIOL

Code Z19316
TO USE WITH THE KIT CODE GC19310

PURPOSE

These calibrators in urine are used to calibrate the GC system for the quantitative measurement of analytes contained therein. Each analyte has a concentration that allows the construction of a curve on a point. The calibrators are freeze-dried human matrix and should be handled as if they were a real patient sample.

RECONSTITUTION


Remove the metal seal and rubber stopper from the vial. Add exactly 10 ml of HPLC grade H₂O in the vial. Replace the rubber stopper, shake and let stand for 5 to 10 minutes. Before use, mix by inverting the vial to dissolve the material until a homogeneous clear solution.

STORAGE AND STABILITY

The calibrators are stable for 36 months from the date of preparation if stored at 2-8 ° C. After reconstitution are stable 7 days at 2-8 ° C and 3 months at -20 ° C. Do not use after the expiry date.

PRECAUTIONS

These calibrators in human matrix should be treated with care and treated as potentially infectious.

LOT		

PACKAGINGS AVAILABLE:




- Z19316 CALIBRATOR IN URINE FOR trans 1,2 CYCLOHEXANEDIOL

4 x 10 ml

N° 003

07/2024

CONCENTRATIONS**CODE****Z19316**

<i>ANALYTE</i>	<i>MEASUREMENT UNIT</i>	<i>VALUE</i>
trans 1,2-CYCLOHEXANEDIOL	mg/l	

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.

**FOR *IN VITRO* DIAGNOSTIC USE ONLY**

EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com