



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS  
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

## CONTROL LYOPHIL. IN URINE MULTIPARAMETRIC FOR ABUSE DRUGS – LEVELS 1 and 2

*(3,4-MDMA, 3,4-MDA, 3,4-MDE, Amphetamine, Methamphetamine, Ketamine, MBDB,  
Ephedrine, Pseudoephedrine, Norpseudoephedrine, Beg, Cocaethylene, Methadone)*

**Code CC43019**

### OBJECTIVE

These controls in urine are used for internal quality control and serve to keep under control the accuracy and precision of the analytical procedures dedicated to the quantitative determination of the analytes contained therein. These controls are lyophilized in human matrix and are available in two different concentration range.

### RECONSTITUTION




Remove the metal seal and rubber stopper from the vial. Add exactly 5 ml of HPLC grade H<sub>2</sub>O in the vial. Replace the rubber stopper, shake and let stand for 5 to 10 minutes. Before use, mix by inverting the vial to dissolve the material until a homogeneous clear solution.

### STORAGE AND STABILITY

The controls are stable for 36 months from the date of preparation if stored at 2-8 ° C. After reconstitution are stable 7 days at 2-8 ° C and 1 month at -20 ° C. Do not use after the expiry date.

### PRECAUTIONS

These controls in human matrix should be treated with care and treated as potentially infectious.

		
019	March 2024	March 2027




### PACKAGING AVAILABLE:

- CC43019 CONTROL IN URINE FOR DRUGS OF ABUSE - LEVELS 1 and 2    2 x 6 x 5 ml

Release N° 005

July 2019

CONCENTRATIONS  
CODE CC43019

		
019	March 2024	March 2027

<i>ANALYTE</i>	<i>MEASURE MENT UNIT</i>	<i>L1</i>	<i>RANGE L1</i>	<i>L2</i>	<i>RANGE L2</i>
3,4 – MDMA	ng/ml	120,0	96,0 – 144,0	660,0	561,0 – 759,0
3,4 – MDA	ng/ml	130,0	104,0 – 156,0	580,0	493,0 – 667,0
3,4 – MDE	ng/ml	120,0	96,0 – 144,0	500,0	425,0 – 575,0
AMPHETAMINE	ng/ml	140,0	112,0 – 168,0	600,0	510,0 – 690,0
METHAMPHETAMINE	ng/ml	140,0	112,0 – 168,0	500,0	425,0 – 575,0
BEG	ng/ml	40,0	32,0 – 48,0	240,0	204,0 – 276,0
COCAETHYLENE	ng/ml	50,0	40,0 – 60,0	340,0	289,0 – 391,0
METHADONE	ng/ml	100,0	80,0 – 120,0	520,0	442,0 – 598,0
MBDB	ng/ml	140,0	112,0 – 168,0	700,0	595,0 – 805,0
KETAMINE	ng/ml	100,0	80,0 – 120,0	650,0	552,5 – 747,5
EPHEDRINE	ng/ml	100,0	80,0 – 120,0	680,0	578,0 – 782,0
PSEUDOEPHEDRINE	ng/ml	140,0	112,0 – 168,0	800,0	680,0 – 920,0
NORPSEUDOEPHEDRINE	ng/ml	140,0	112,0 – 168,0	600,0	510,0 – 690,0

\* The procedures for production and validation / quality control of the individual lot shall be performed using the chromatographic method LC / MS / MS validated according to international guidelines of the Food and Drug Administration (FDA-May 2001) and the International Conference of Harmonisation of technical requirements for registration of pharmaceuticals for human use (ICH-Q2-R1-Nov 2005). The mean value and standard deviation were obtained in the validation phase of batch processing for level four bottles in duplicate.

This product fulfils all the requirements of Directive 98/79/EC of 27/10/1998 on *in vitro* diagnostic medical devices (IVD). The declaration of conformity is available upon request.



## FOR *IN VITRO* DIAGNOSTIC USE ONLY

produced by

COMPANY WITH  
MANAGEMENT SYSTEM  
CERTIFIED BY DNV  
= ISO 9001 =  
= ISO 13485 =



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS  
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

EUREKA srl – LAB DIVISION ,  
Via Enrico Fermi 25, 60033 Chiaravalle (AN) ITALY  
Tel. + 39 071 74 50 790  
Fax. + 39 071 74 96 579  
[info@eurekaone.com](mailto:info@eurekaone.com)  
[www.eurekakit.com](http://www.eurekakit.com)

