

CALIBRATOR LYOPHIL. IN PLASMA FOR HOMOCYSTEINE

Code LC09016
TO USE WITH THE KIT CODE LC09010-15

OBJECTIVE

The calibrators in plasma are used to calibrate the HPLC system for the quantitative measurement of analytes contained therein. Each analyte has a concentration that allows the construction of a calibration curve of five points. The calibrators are freeze-dried human matrix and should be handled as if they were a real patient sample.

RECONSTITUTION

Remove the metal seal and rubber stopper from the vial. Add exactly 1 ml of HPLC grade H₂O in the vial. Replace the rubber stopper, shake and let stand for 5 to 10 minutes. Before use, mix by inverting the vial to dissolve the material until a homogeneous clear solution.

STORAGE AND STABILITY

The calibrators are stable for 36 months from the date of preparation if stored at -20° C. After reconstitution are stable 7 days at 2-8 ° C and 6 months at -20 ° C. Do not use after the expiry date.

PRECAUTIONS

These controls in human matrix should be treated with care and treated as potentially infectious.

		
006	Oct 2020	Oct 2023

PACKAGINGS AVAILABLE:

- LC09016 CALIBRATOR IN PLASMA FOR HOMOCYSTEINE 5 x 2 x 1 ml

CONCENTRATIONS**CODE****LC09016**

		
006	Oct 2020	Oct 2023

<i>ANALYTE</i>	<i>MEASUREMENT UNIT</i>	<i>VALUE</i>
HOMOCYSTEINE - 1 -	μmoli/l	0,7
HOMOCYSTEINE - 2 -	μmoli/l	4,7
HOMOCYSTEINE - 3 -	μmoli/l	8,8
HOMOCYSTEINE - 4 -	μmoli/l	17,4
HOMOCYSTEINE - 5 -	μmoli/l	42,5

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 and D.L. n.332 of 08/09/2000 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.

**FOR IN VITRO DIAGNOSTIC USE ONLY**

produced by



The ONLY Italian Company that Researches, Develops and Produces
diagnostic KITS for HPLC, GC, GC-MS, LC-MS/MS

COMPANY WITH
MANAGEMENT SYSTEM
CERTIFIED BY DNV
= ISO 9001 =
= ISO 13485 =

EUREKA srl – LAB DIVISION ,
Via Enrico Fermi 25, 60033 Chiaravalle (AN) ITALY
Tel. + 39 071 74 50 790 Fax. + 39 071 74 96 579
info@eurekaone.com, info@eurekakit.com
www.eurekaone.com, www.eurekakit.com

Release N° 003

September 2019