



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CALIBRATOR LYOPHIL. IN PLASMA FOR ANTIARRHYTHMICS
(Amiodarone, Desethylamiodarone, Flecainide, Propafenone, Sotalol, Quinidine)

Code LC99116

OBJECTIVE

These plasma calibrators are used to calibrate the LC system in the quantitative determination of the analytes they contain. Each analyte has a concentration that allows the construction of a seven-point calibration curve. These freeze-dried calibrators are in human matrix and should be handled as if they were a real patient sample.

RECONSTITUTION

Remove the metal seal and rubber cap from the vial. Add exactly 1 ml of HPLC grade water to the vial. Replace the rubber cap, shake and leave to stand for 5 - 10 minutes. Before use shake the vial by inversion to dissolve the material until a clear and homogeneous solution is obtained.

STORAGE AND STABILITY

The calibrators are stable for 36 months from the date of preparation if stored at 2-8 ° C. After reconstitution are stable 7 days at 2-8 ° C and 1 month at -20 ° C. Do not use after the expiry date.

PRECAUTIONS

These calibrators in human matrix should be treated with care and treated as potentially infectious.

LOT		

AVAILABLE PACKAGING:

- LC99116 CALIBRATOR IN PLASMA FOR ANTIARRHYTHMICS 7 x 2 x 1 ml

N° 002	10/2024
--------	---------

CONCENTRATIONS
CODE LC99116

LOT		

ANALYTE	MEASUREMENT UNIT	C0	C1	C2	C3	C4	C5	C6
AMIODARONE	mg/l							
DESETHYLAMIODARONE	mg/l							
FLECAINIDE	mg/l							
PROPAFENONE	mg/l							
SOTALOL	mg/l							
QUINIDINE	mg/l							

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR IN VITRO DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com