



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CONTROL LYOPHIL. IN URINE FOR T,T-MUCONIC ACID and S-PHENYLMERCAPTURIC ACID - LEVELS 1 and 2

Code LC23019

OBJECTIVE

These urine controls are used for internal quality assurance in clinical-chemical laboratories and are made to keep under control the care and precision of the analytical procedure in charge to the quantitative determination of in the analytes in them included. These lyophilised controls are based on human urine and are available with mean values in the normal as well as in the pathological range. The controls have to be prepared like patient samples in one series of analyses.

RECONSTITUTION

Remove the metal seal and rubber stopper from the vial. Add exactly 1 ml of HPLC grade H₂O in the vial. Replace the rubber stopper, shake and let stand for 5 to 10 minutes. Before use, mix by inverting the vial to dissolve the material until a homogeneous clear solution.

STORAGE AND STABILITY

The calibrators are stable for 36 months from the date of preparation if stored at -20 ° C. After reconstitution are stable 7 days at 2-8 ° C and 6 months at -20 ° C. Do not use after the expiry date.

PRECAUTIONS

These controls in human matrix should be treated with care and treated as potentially infectious.

		

AVAILABLE PACKAGINGS:

- LC23019 CONTROL IN URINE FOR T,T-MUCONIC/S-PHENYLMERCAPTURIC AC. - LEVELS 1 and 2
2 x 5 x 1 ml

N° 004	06/2024
--------	---------

CONCENTRATIONS
CODE LC23019

LOT		

<i>ANALYTE</i>	<i>MEASUREMENT UNIT</i>	<i>MEDIUM VALUE L1</i>	<i>RANGE L1</i>	<i>MEDIUM VALUE L2</i>	<i>RANGE L2</i>
T,T-MUCONIC ACID	µg/l				
S-PHENYLMERCAPTURIC ACID	µg/l				

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR IN VITRO DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com