



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS  
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

**CONTROL LYOPHIL. IN SERUM FOR BENZODIAZEPINES – LEVELS 1, 2 and 3**  
(*Bromazepam, Chlordemethyldiazepam, Chlordiazepoxide, Clobazam, Clonazepam, Demoxepam, Diazepam, Flunitrazepam, Flurazepam, Lorazepam, Midazolam, Nitrazepam, Norclobazam, Nordiazepam, Oxazepam, Trazodone, Triazolam, Zopiclone*)

**Code LC80019**

**PURPOSE**

These plasma controls are used for internal quality assurance in clinical-chemical laboratories and are made to keep under control the care and precision of the analytical procedure in charge to the quantitative determination of in the analytes in them included. These lyophilised controls are based on human plasma and are available with mean values in the normal as well as in the pathological range. The controls have to be prepared like patient samples in one series of analyses.

**RECONSTITUTION**

Remove the metal seal and rubber stopper from the vial. Add exactly 1 ml of HPLC grade H<sub>2</sub>O in the vial. Replace the rubber stopper, shake and let stand for 5 to 10 minutes. Before use, mix by inverting the vial to dissolve the material until a homogeneous clear solution.

**STORAGE AND STABILITY**

The calibrators are stable for 36 months from the date of preparation if stored at 2-8 ° C. After reconstitution are stable 7 days at 2-8 ° C and 1 month at -20 ° C. Do not use after the expiry date.

**PRECAUTIONS**

These controls in human matrix should be treated with care and treated as potentially infectious.

|                                                                                     |                                                                                     |                                                                                     |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |
|                                                                                     |                                                                                     |                                                                                     |

**PACKAGING AVAILABLE:**

- LC80019 CONTROL IN SERUM FOR BENZODIAZEPINES - LEVELS 1, 2 and 3      3 x 5 x 1 ml

|        |         |
|--------|---------|
| N° 003 | 08/2024 |
|--------|---------|

CONCENTRATIONS  
CODE LC80019

|     |                                                                                     |                                                                                     |
|-----|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| LOT |  |  |
|     |                                                                                     |                                                                                     |

| ANALYTE               | MEASUREMENT UNIT | MEDIUM VALUE L1 | RANGE L1 | MEDIUM VALUE L2 | RANGE L2 | MEDIUM VALUE L3 | RANGE L3 |
|-----------------------|------------------|-----------------|----------|-----------------|----------|-----------------|----------|
| BROMAZEPAM            | µg/l             |                 |          |                 |          |                 |          |
| CHLORDEMETHYLDIAZEPAM | µg/l             |                 |          |                 |          |                 |          |
| CHLORDIAZEPOXIDE      | µg/l             |                 |          |                 |          |                 |          |
| CLOBAZAM              | µg/l             |                 |          |                 |          |                 |          |
| CLONAZEPAM            | µg/l             |                 |          |                 |          |                 |          |
| DEMOXEPAM             | µg/l             |                 |          |                 |          |                 |          |
| DIAZEPAM              | µg/l             |                 |          |                 |          |                 |          |
| FLUNITRAZEPAM         | µg/l             |                 |          |                 |          |                 |          |
| FLURAZEPAM            | µg/l             |                 |          |                 |          |                 |          |
| LORAZEPAM             | µg/l             |                 |          |                 |          |                 |          |
| MIDAZOLAM             | µg/l             |                 |          |                 |          |                 |          |
| NITRAZEPAM            | µg/l             |                 |          |                 |          |                 |          |
| NORCLOBAZAM           | µg/l             |                 |          |                 |          |                 |          |
| NORDIAZEPAM           | µg/l             |                 |          |                 |          |                 |          |
| OXAZEPAM              | µg/l             |                 |          |                 |          |                 |          |
| TRAZODONE             | µg/l             |                 |          |                 |          |                 |          |
| TRIAZOLAM             | µg/l             |                 |          |                 |          |                 |          |
| ZOPICLONE             | µg/l             |                 |          |                 |          |                 |          |

This product fulfils all the requirements of Directive 98/79/EC of 27/10/1998 on *in vitro* diagnostic medical devices (IVD). The declaration of conformity is available upon request.



**FOR *IN VITRO* DIAGNOSTIC USE ONLY**



EUREKA S.R.L. LAB DIVISION  
Via M. D'Antona 28  
60033 Chiaravalle (AN) - Italy  
Tel +39 071 7450790  
[eureka-support@sentinel.it](mailto:eureka-support@sentinel.it)  
[www.eurekakit.com](http://www.eurekakit.com)