



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CALIBRATOR LYOPHIL. IN PLASMA FOR ANTITUBERCOLOSIS DRUGS
(Isoniazid, Rifampicin, Ethambutol)

Code LC79116
TO USE WITH KIT CODE LC79110

PURPOSE

These plasma calibrators are used to calibrate the LC system for the quantitative measurement of analytes contained therein. Each analyte has a concentration that allows the construction of a curve on 6 points. The calibrators are freeze-dried human matrix and should be handled as if they were a real patient sample.

RECONSTITUTION

Remove the metal seal and rubber stopper from the vial. Add exactly 1 ml of HPLC grade H₂O in the vial. Replace the rubber stopper, shake and let stand for 5 to 10 minutes. Before use, mix by inverting the vial to dissolve the material until a homogeneous clear solution.

STORAGE AND STABILITY

The calibrators are stable for 36 months from the date of preparation if stored at 2-8 ° C. After reconstitution are stable 7 days at 2-8 ° C and 2 months at -20 ° C. Do not use after the expiry date.

PRECAUTIONS

These calibrators in human matrix should be treated with care and treated as potentially infectious.

LOT		

PACKAGING AVAILABLE:

- LC79116 CALIBRATOR IN PLASMA FOR ANTITUBERCOLOSIS DRUGS 6 x 2 x 1 ml

N° 002	08/2024
--------	---------

CONCENTRATIONS
CODE LC79116

LOT		

ANALYTE	MEASUREMENT UNIT	C0	C1	C2	C3	C4	C5
ETHAMBUTOL	mg/l						
ISONIAZID	mg/l						
RIFAMPICIN	mg/l						

This product fulfils all the requirements of Directive 98/79/EC of 27/10/1998 on *in vitro* diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR *IN VITRO* DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com