



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CALIBRATOR FOR PHYTOCANNABINOIDS

(Tetrahydrocannabinol [THC], Cannabidiol [CBD], Cannabinol [CBN], Cannabigerol [CBG],
Tetrahydrocannabinolic Acid [THCA], Cannabidiolic Acid [CBDA], Tetrahydrocannabivarin [THCV])

Code GC88816

OBJECTIVE



These calibrators are used to calibrate the GC system in the quantitative determination of the analytes they contain. Each analyte has a concentration such that a six-point calibration curve can be constructed.

STORAGE AND STABILITY

Once the vial is opened, take the necessary volume, close it with the unused reagent and store at -20°C. Stable 6 months. Do not use after expiry date.

PRECAUTIONS

These calibrators should be treated with care and considered as potentially infectious.

LOT		



AVAILABLE PACKAGINGS:

- GC88816 CALIBRATOR FOR PHYTOCANNABINOIDS 6 x 2 x 600 µl

N° 005

09/2025

CONCENTRATIONS
CODE GC88816

LOT		

ANALYTE	MEASUREMENT UNIT	C0	C1	C2	C3	C4	C5
THC	mg/ml						
CBD	mg/ml						
CBN	mg/ml						
THCA	mg/ml						
THCV	mg/ml						
CBG	mg/ml						
CBDA	mg/ml						

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR IN VITRO DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com