



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CONTROLS IN SALIVA FOR CORTISOL AND CORTISONE – LEVELS 1, 2 and 3

Code LC91119

OBJECTIVE

These controls in saliva are used for internal quality control and serve to keep under control the accuracy and precision of the analytical procedures dedicated to the quantitative determination of the analytes contained therein.

STORAGE AND STABILITY

Once you open the vial, withdrawn the volume needed, recap and stored at  -15°C to -25°C . Stable 6 months. Do not use after the expiration date.

PRECAUTIONS

These controls should be treated with care and treated as potentially infectious.

LOT		

AVAILABLE PACKAGING:

- LC91119 CONTROLS IN SALIVA FOR CORTISOL AND CORTISONE – LEVELS 1, 2 and 3
3 x 5 x 3 ml

N° 003	06/2024
--------	---------

CONCENTRATIONS
CODE LC91119

LOT		

ANALYTE	MEASUREMENT UNIT	L1	RANGE	L2	RANGE	L3	RANGE
CORTISOL	ng/ml						
CORTISONE	ng/ml						

* The procedures for production and validation / quality control of every batch shall be performed using the chromatographic method LC / MS / MS validated according to international guidelines of the Food and Drug Administration (FDA-May 2001) and the International Conference of Harmonisation of technical requirements for registration of pharmaceuticals for human use (ICH-Q2-R1-Nov 2005). The mean value and standard deviation were obtained during the validation phase of batch processing for level four bottles in duplicate.

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR IN VITRO DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com