



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CALIBRATOR FOR PHYTOCANNABINOIDS

*(Tetrahydrocannabinol [THC], Cannabidiol [CBD], Cannabinol [CBN], Cannabigerol [CBG],
Tetrahydrocannabinolic Acid [THCA], Cannabidiolic Acid [CBDA], Tetrahydrocannabivarin [THCV])*

Code LC88816

OBJECTIVE

These calibrators are used to calibrate the HPLC system in the quantitative determination of the analytes they contain. Each analyte has a concentration such that a six-point calibration curve can be constructed.

STORAGE AND STABILITY

Once the vial is opened, take the necessary volume, close it with the unused reagent and store at -20°C. Stable 6 months. Do not use after expiry date.

PRECAUTIONS

These calibrators should be treated with care and considered as potentially infectious.

LOT		

AVAILABLE PACKAGINGS:

- LC88816 CALIBRATOR FOR PHYTOCANNABINOIDS 6 x 2 x 300 µl

N° 004	07/2024
--------	---------

CONCENTRATIONS
CODE LC88816

		

<i>ANALYTE</i>	<i>MEASUREMENT UNIT</i>	<i>C0</i>	<i>C1</i>	<i>C2</i>	<i>C3</i>	<i>C4</i>	<i>C5</i>
THC CBD CBN THCA THCV CBG CBDA	mg/ml						

This product fulfills all the requirements of Directive 98/79/EC of 27/10/1998 on in vitro diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR *IN VITRO* DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
 Via M. D'Antona 28
 60033 Chiaravalle (AN) - Italy
 Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com