



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

Mobile Phase M for Homocysteine in plasma

Code Z09611

This Mobile Phase can be used as indicated in IFU kit code Z09610-Z09710
as Reagent M (Mobile Phase M)

Packaging **4 x 500 ml**

The Mobile Phase is stable for 3 years from the preparation date at 2-8 °C.

This product fulfils all the requirements of Directive 98/79/EC of 27/10/1998 on *in vitro* diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR *IN VITRO* DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com