



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

CALIBRATOR FOR VOC IN URINE

(1-BUTANOL, ACETONE, BENZENE, CYCLOHEXANE, CYCLOHEXANOL, CYCLOHEXANONE, CHLOROFORM, DICHLOROMETHANE, ETHANOL, ISOBUTANOL, ISOPROPANOL, METHANOL, METHYL ETHYL KETONE, METHYL ISOBUTHYL KETONE, PERCHLOROETHYLENE, STYRENE, TOLUENE, o-m-p-XYLENE)

Code GC77016
TO USE WITH KIT CODE GC77010

OBJECTIVE



These calibrators are used to calibrate the GC/MS system in the quantitative determination of the analyte contained in them. Each analyte has a concentration such that a 7-point calibration curve can be constructed. These calibrators are in human matrix and should be handled as if they were a real patient sample.

STORAGE AND STABILITY

The calibrators are stable 36 months from the date of preparation if stored at -25/-15°C. Do not use them after the expiry date.

PRECAUTIONS

These calibrators in human matrix must be treated with care and considered as potentially infectious.

LOT		

AVAILABLE PACKAGINGS:



GC77016 CALIBRATOR FOR VOC IN URINE

7 x 2 x 10 ml

N° 005

11/2025

CONCENTRATIONS
CODE GC77016

LOT		

ANALYTE	MEASUREMENT UNIT	C0	C1	C2	C3	C4	C5	C6
1-BUTANOL	mg/l							
ACETONE	mg/l							
BENZENE	mg/l							
CYCLOHEXANE	mg/l							
CYCLOHEXANOL	mg/l							
CTCLOHEXANONE	mg/l							
CHLOROFORM	mg/l							
DICHLOROMETHANE	mg/l							
ETHANOL	mg/l							
ISOBUTANOL	mg/l							
ISOPROPANOL	mg/l							
METHANOL	mg/l							
MEK	mg/l							
MIBK	mg/l							
PERCHLOROETHYLENE	mg/l							
STYRENE	mg/l							
TOLUENE	mg/l							
o-XYLENE	mg/l							
∑ m-p-XYLENE	mg/l							

This product fulfils all the requirements of Directive 98/79/EC of 27/10/1998 on *in vitro* diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR *IN VITRO* DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com