



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS  
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

## CONTROL FOR VOC IN URINE

### Levels 1 and 2

(1-BUTANOL, ACETONE, BENZENE, CYCLOHEXANE, CYCLOHEXANOL, CYCLOHEXANONE, CHLOROFORM, DICHLOROMETHANE, ETHANOL, ISOBUTANOL, ISOPROPANOL, METHANOL, METHYL ETHYL KETONE, METHYL ISOBUTHYL KETONE, PERCHLOROETHYLENE, STYRENE, TOLUENE, o-m-p-XYLENE)

### Code GC77019

#### OBJECTIVE



These controls are used for internal quality control and serve to control the accuracy and precision of analytical procedures dedicated to the quantitative determination of the analytes they contain. These controls are in human matrix and are available in two different concentration ranges.

#### STORAGE AND STABILITY

The controls are stable 36 months from the date of preparation if stored at -25/-15°C. Do not use them after the expiry date.

#### PRECAUTIONS

These controls in human matrix must be treated with care and considered as potentially infectious.



LOT		

#### AVAILABLE PACKAGING:

- GC77019 CONTROL FOR VOC IN URINE - Levels 1 and 2

2 x 10 x 10 ml

CONCENTRATIONS  
CODE GC77019

LOT		

<b>ANALYTE</b>	<b>MEASUREMENT UNIT</b>	<b>MEDIUM VALUE L1</b>	<b>RANGE L1</b>	<b>MEDIUM VALUE L2</b>	<b>RANGE L2</b>
1-BUTANOL	mg/l				
ACETONE	mg/l				
BENZENE	mg/l				
CYCLOHEXANE	mg/l				
CYCLOHEXANOL	mg/l				
CYCLOHEXANONE	mg/l				
CHLOROFORM	mg/l				
DICHLOROMETHANE	mg/l				
ETHANOL	mg/l				
ISOBUTANOL	mg/l				
ISOPROPANOL	mg/l				
METHANOL	mg/l				
MEK	mg/l				
MIBK	mg/l				
PERCHLOROETHYLENE	mg/l				
STYRENE	mg/l				
TOLUENE	mg/l				
o-XYLENE	mg/l				
Σ m-p-XYLENE	mg/l				

This product fulfils all the requirements of Directive 98/79/EC of 27/10/1998 on *in vitro* diagnostic medical devices (IVD). The declaration of conformity is available upon request.



**FOR *IN VITRO* DIAGNOSTIC USE ONLY**



EUREKA S.R.L. LAB DIVISION  
Via M. D'Antona 28  
60033 Chiaravalle (AN) - Italy  
Tel +39 071 7450790  
[eureka-support@sentinel.it](mailto:eureka-support@sentinel.it)  
[www.eurekakit.com](http://www.eurekakit.com)