



AN ITALIAN COMPANY THAT RESEARCHES, DEVELOPS AND PRODUCES READY TO USE KITS
USING CHROMATOGRAPHIC TECHNIQUES FOR CLINICAL LABORATORIES

STARTER KIT FOR DELTA-9-THC AND DELTA 9 THC-COOH IN HAIR BY LC/MS

Code SK53010

OBJECTIVE

The Starter Kit is needed to identify the transitions and collision energies of molecules on a mass detector attached to a chromatographic system. It should be used only for the qualitative analysis of delta-9-THC.

STORE AND STABILITY

The Starter Kit is stable for 12 months from the date of preparation if stored at -25/-15 °C. Do not use after expiry date.

LOT		

AVAILABLE PACKAGING:

- SK53010 STARTER KIT FOR DELTA-9-THC AND DELTA 9 THC-COOH IN HAIR BY LC/MS 4 x 1 ml

This product fulfils all the requirements of Directive 98/79/EC of 27/10/1998 on *in vitro* diagnostic medical devices (IVD). The declaration of conformity is available upon request.



FOR *IN VITRO* DIAGNOSTIC USE ONLY



EUREKA S.R.L. LAB DIVISION
Via M. D'Antona 28
60033 Chiaravalle (AN) - Italy
Tel +39 071 7450790
eureka-support@sentinel.it
www.eurekakit.com

CODE

SK53010

LOT		

Abuse Drugs
Delta-9-THC
Delta-9-THC-COOH
Delta-9-THC-D3
Delta-9-THC-COOH-D9